

Product Datasheet

Histone H4 (Di Methyl Lys20) antibody (orb764083)

Description	Rabbit polyclonal antibody to Histone H4 (Di Methyl Lys20)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide of Histone H4 (Di Methyl Lys20)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14
Expiration Date	12 months from date of receipt.