

Product Datasheet

ASK 1 (phospho-Ser966) antibody (orb764139)

Description Rabbit polyclonal antibody to ASK 1 (phospho-Ser966)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human ASK1 around the phosphorylation site of Ser966. AA range:932-981

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 155


Expiration Date 12 months from date of receipt.



Western Blot analysis of various cells u...



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from 29...