



Product Datasheet

Chk2 (phospho-Thr68) antibody (orb764163)


Description	Rabbit polyclonal antibody to Chk2 (phospho-Thr68)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Chk2 around the phosphorylation site of Thr68. AA range:35-84
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Hela cell...



Immunofluorescence analysis of rat-heart...



Immunofluorescence analysis of rat-lung ...