


Product Datasheet


Cyclin E1 (phospho-Thr395) antibody (orb764170)

Description	Rabbit polyclonal antibody to Cyclin E1 (phospho-Thr395)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Cyclin E1 around the phosphorylation site of Thr395. AA range:361-410
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48
Expiration Date	12 months from date of receipt.




Western Blot analysis of various cells

U...



Western Blot analysis of HeLa cells using...



Enzyme-Linked Immunosorbent Assay (Phosp...