

Product Datasheet

Flk-1 (phospho-Tyr1214) antibody (orb764189)

Description Rabbit polyclonal antibody to Flk-1 (phospho-Tyr1214)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human VEGFR2 around the phosphorylation site of Tyr1214. AA range:1180-1229

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 152


Expiration Date 12 months from date of receipt.



Western Blot analysis of various cells u...



Enzyme-Linked Immunosorbent Assay (Phosp...)



Immunohistochemistry analysis of paraffi...