



Product Datasheet

Histone H2A.X (phospho-Ser139) antibody (orb764201)


Description	Rabbit polyclonal antibody to Histone H2A.X (phospho-Ser139)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Histone H2A.X around the phosphorylation site of Ser139. AA range:94-143
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15 19
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from 1)...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...