



Product Datasheet

Integrin beta1 (phospho-Thr788) antibody (orb764213)


Description	Rabbit polyclonal antibody to Integrin beta1 (phospho-Thr788)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Integrin beta1 around the phosphorylation site of Thr788. AA range:749-798
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	135
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from 1)...



Western Blot analysis of various cells u...



Enzyme-Linked Immunosorbent Assay (Phosp...