


Product Datasheet


Stat6 (phospho-Tyr641) antibody (orb764284)

Description	Rabbit polyclonal antibody to Stat6 (phospho-Tyr641)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human STAT6 around the phosphorylation site of Tyr641. AA range:608-657
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	94
Expiration Date	12 months from date of receipt.




Western blot analysis of Stat6 (phospho-Tyr641)

...



Western Blot analysis of various cells

u...



Enzyme-Linked Immunosorbent Assay (Phospho-Tyr641)