



Product Datasheet

PKD2 (phospho-Ser876) antibody (orb764307)


Description	Rabbit polyclonal antibody to PKD2 (phospho-Ser876)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PKD2 around the phosphorylation site of Ser876. AA range:829-878
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	96
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells u...



Western Blot analysis of HELA cells usin...



Enzyme-Linked Immunosorbent Assay (Phosp...