



Product Datasheet

RPA p32 (phospho-Ser33) antibody (orb764314)


Description	Rabbit polyclonal antibody to RPA p32 (phospho-Ser33)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RFA2 around the phosphorylation site of Ser33. AA range:1-50
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
various
cells u...



Western
Blot
analysis of
MCF-7 cells
usi...



Western
blot
analysis of
lysates from
NI...