



Product Datasheet

EphB1/2 (phospho-Tyr594/604) antibody (orb764324)


Description	Rabbit polyclonal antibody to EphB1/2 (phospho-Tyr594/604)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human EPHB1/2 around the phosphorylation site of Tyr594/604. AA range:561-610
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells u...



Western Blot analysis of Hela cells usin...



Immunofluorescence analysis of HUVEC cel...