



Product Datasheet

NF kappa B-p105/p50 (phospho-Ser337) antibody (orb764341)


Description	Rabbit polyclonal antibody to NF kappa B-p105/p50 (phospho-Ser337)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized phospho-peptide around the phosphorylation site of human NFkB-p105/p50 (phospho Ser337)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105+50
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells u...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...