



Product Datasheet

ADAMTS-1 antibody (orb764474)

Description	Rabbit polyclonal antibody to ADAMTS-1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from ADAMTS-1 . at AA range: 160-240
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells
U...



Western Blot analysis of HUVEC cells
usi...