



Product Datasheet

ALDH1B1 antibody (orb764512)


Description	Rabbit polyclonal antibody to ALDH1B1
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human ALDH1B1. AA range:311-360
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells
u...



Western Blot analysis of COS7 cells
usin...



Western blot analysis of lysates from He...