


Product Datasheet


AMPKalpha1 antibody (orb764523)

Description	Rabbit polyclonal antibody to AMPKalpha1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human AMPK1. AA range:451-500
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of rat-lung

...



Immunofluorescence analysis of rat-lung

...



Western Blot analysis of mouse-lung cell...