



Product Datasheet

ANT3 antibody (orb764535)


Description	Rabbit polyclonal antibody to ANT3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC25A6. AA range:121-170
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells
U...



Western blot analysis of lysates from HepG2 and HeLa cells



Western blot analysis of the lysates from HepG2 cells