



Product Datasheet

APRIL antibody (orb764560)


Description	Rabbit polyclonal antibody to APRIL
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human TALL-2. AA range:151-200
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23
Expiration Date	12 months from date of receipt.



Western Blot
analysis of various
cells u...



Western Blot
analysis of COLO205
cells u...



Immunofluorescence
analysis of HeLa
cell...