



Product Datasheet

ARK-2 antibody (orb764576)


Description	Rabbit polyclonal antibody to ARK-2
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human AurB. AA range:198-247
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39
Expiration Date	12 months from date of receipt.



Western Blot
analysis of various
cells u...



Immunofluorescence
analysis of HeLa
cell...



Western blot
analysis of lysates
from CO...