


Product Datasheet

C6 antibody (orb764685)


Description	Rabbit polyclonal antibody to C6
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human C6. AA range:221-270
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells u...



Western blot analysis of C6 Antibody. Th...



Immunohistochemistry analysis of paraff...