


Product Datasheet

CD44 antibody (orb7647)

Description	Mouse monoclonal antibody to CD44
Species/Host	Mouse
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	Stimulated human leukocytes
Target	CD44
Preservatives	Purified antibody in PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	p16070
Entrez	960, 101008690
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try at: Flow Cyt 1:10, ICC/IF 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
using CD44
a...