



Product Datasheet

Chk2 antibody (orb764849)


Description	Rabbit polyclonal antibody to Chk2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CHEK2. AA range:35-84
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61
Expiration Date	12 months from date of receipt.



Western blot analysis of Chk2 Polyclonal...



Western Blot analysis of various cells u...



Western Blot analysis of HuvEc cells usi...