



Product Datasheet

Collagen II antibody (orb764897)

Description	Rabbit polyclonal antibody to Collagen II
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Collagen II. AA range:101-150
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140
Expiration Date	12 months from date of receipt.




Western blot analysis of lysates from 1)...



Immunofluorescence analysis of rat-lung

...



Immunofluorescence analysis of rat-lung

...