



Product Datasheet

CYP2A7 antibody (orb764973)


Description	Rabbit polyclonal antibody to CYP2A7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from the C-terminal region of human CYP2A7.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells u...



Western Blot analysis of 22RV1 cells usi...



Immunohistochemical analysis of paraffin...