

## Product Datasheet

### MST1R (phospho-Ser1349) antibody (FITC) (orb7650)

|                            |                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to MST1R (phospho-Ser1349) antibody conjugated to FITC                                      |
| <b>Species/Host</b>        | Rabbit                                                                                                        |
| <b>Conjugation</b>         | FITC                                                                                                          |
| <b>Tested Applications</b> | FC, IF                                                                                                        |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from mouse MST1R around the phosphorylation site of Ser1349 |
| <b>Target</b>              | phospho-Ron (Ser1349)                                                                                         |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                                         |
| <b>Form/Appearance</b>     | Lyophilized or Liquid                                                                                         |
| <b>Concentration</b>       | - 100 µg (in 100 µl): 1 mg/ml                                                                                 |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                              |
| <b>Note</b>                | For research use only                                                                                         |
| <b>Isotype</b>             | IgG                                                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                                                    |
| <b>MW</b>                  | 30/117/148kDa                                                                                                 |
| <b>Uniprot ID</b>          | <a href="#">Q04912</a>                                                                                        |
| <b>Dilution Range</b>      | Flow-Cyt=3ug/Test,IF=1:100-500                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                               |