



Product Datasheet

Elf-5 antibody (orb765126)


Description	Rabbit polyclonal antibody to Elf-5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human ELF5. AA range:191-240
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells U...



Western blot analysis of lysates from Ju...



Western blot analysis of the lysates fro...