


Product Datasheet


Flk-1 antibody (orb765226)

Description	Rabbit polyclonal antibody to Flk-1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human VEGFR2. AA range:917-966
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of rat-lung

...



Immunofluorescence analysis of rat-lung

...



Immunofluorescence analysis of rat-kidney...