



Product Datasheet

GLP-1 antibody (orb765307)


Description	Rabbit polyclonal antibody to GLP-1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human GLP1R. AA range:101-150
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using GLP-1 antibody (orb765307).



Western Blot analysis of various cells using GLP-1 antibody (orb765307).



Immunofluorescence analysis of HeLa cells using GLP-1 antibody (orb765307).