



Product Datasheet

GNL1 antibody (orb765316)


Description	Rabbit polyclonal antibody to GNL1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human GNL1. AA range:61-110
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	69
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells
u...



Western Blot analysis of K562 cells
usin...



Western blot analysis of lysates from
Ju...