



Product Datasheet

IL17RC antibody (orb765484)


Description	Rabbit polyclonal antibody to IL17RC
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human IL17RC. AA range:721-770
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70
Expiration Date	12 months from date of receipt.



Western Blot analysis of HuvEc cells using IL17RC antibody



Western Blot analysis of various cells using IL17RC antibody



Western blot analysis of the lysates from HuvEc cells using IL17RC antibody