


Product Datasheet

Niban (Phospho-Ser602) Antibody (orb76555)

Description	Niban (Phospho-Ser602) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BZQ8
Expiration Date	12 months from date of receipt.



Western blot analysis of Niban antibody ...