



Product Datasheet

KIR6.2 antibody (orb765551)


Description	Rabbit polyclonal antibody to KIR6.2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Kir6.2. AA range:190-239
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells
u...



Western Blot analysis of NIH-3T3 cells
u...



Western Blot analysis of various cells
u...