


Product Datasheet

S6 Ribosomal Protein (Phospho-Ser235/236) Antibody (orb76557)

Description	S6 Ribosomal Protein (Phospho-Ser235/236) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P62753
Expiration Date	12 months from date of receipt.



Western blot analysis of S6 Ribosomal Protein.