



Product Datasheet

LATS1/2 antibody (orb765579)


Description	Rabbit polyclonal antibody to LATS1/2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human LATS1/2. AA range:1041-1090
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130-140
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549. 1, ...



Western Blot analysis of mouse-muscle cel...



Immunohistochemistry analysis of paraffi...