



Product Datasheet

MMP-1 antibody (orb765677)


Description	Rabbit polyclonal antibody to MMP-1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MMP-1. AA range:411-460
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54
Expiration Date	12 months from date of receipt.



Western blot analysis of HELA 293T 3T3 I...



Western Blot analysis of HepG2 cells usi...



Western blot analysis of lysates from 1)...