



Product Datasheet

Myomesin-2 antibody (orb765754)


Description	Rabbit polyclonal antibody to Myomesin-2
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MYOM2. AA range:612-661
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	165
Expiration Date	12 months from date of receipt.



Western Blot analysis of HT29 cells using Myomesin-2 antibody



Western Blot analysis of various cells using Myomesin-2 antibody



Western blot analysis of lysates from COS7 and HT29 cells using Myomesin-2 antibody