


Product Datasheet

Acetyl-CoA Carboxylase (Ab-79) Antibody (orb76576)

Description	Acetyl-CoA Carboxylase(Ab-79) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13085
Expiration Date	12 months from date of receipt.



Western blot analysis of Acetyl-CoA Carb...