


## Product Datasheet

### PTH rP antibody (orb76579)


|                            |                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to PTH rP                                                 |
| <b>Species/Host</b>        | Mouse                                                                               |
| <b>Reactivity</b>          | Human                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                        |
| <b>Tested Applications</b> | ELISA, FC, IHC-P, WB                                                                |
| <b>Immunogen</b>           | 23 amino acid peptide (corresponding to amino acids 146-169) of the hPTHrP receptor |
| <b>Target</b>              | Parathyroid Hormone Receptor 1                                                      |
| <b>Preservatives</b>       | Purified antibody containing PBS 0.1% sodium azide                                  |
| <b>Concentration</b>       | 1mg/ml                                                                              |
| <b>Note</b>                | For research use only                                                               |
| <b>Application notes</b>   | P3/X63.Ag8.653 mouse myeloma cells with spleen cells from a BALB/c;                 |
| <b>Isotype</b>             | Ig                                                                                  |
| <b>Clonality</b>           | Monoclonal                                                                          |
| <b>Uniprot ID</b>          | p12272                                                                              |
| <b>Entrez</b>              | <b>5744</b>                                                                         |
| <b>Dilution Range</b>      | Optimal antibody dilution should be determined by titration                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                     |



Immunohistochemical staining of PTHrP an...



Line graph illustrates about the Ag-Ab r...



Western blot analysis of HEK cells using...