


Product Datasheet

Acetylated Lysine antibody (orb76580)

Description	Mouse monoclonal antibody to Acetylated Lysine
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	Acetylated keyhole limpet hemocyanin (KLH) was used as immunogen. The cell line was screened against non acetylated KLH as a control in addition to BSA and histones
Target	Acetyl Lysine
Preservatives	Purified antibody containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try 1:3000-1:10,000 for EIA
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...