


Product Datasheet

ALY antibody (orb76584)

Description	Mouse monoclonal antibody to ALY
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, EM, FC, ICC, IHC-P, IP, WB
Immunogen	Full length recombinant human ALY/Ref fusion protein
Target	ALY/Ref
Preservatives	Purified antibody (from supernatant) containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P22626
Entrez	10189, 21681
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try:IHC-P: Use a concentration of 1 µg/ml.ICC/IF: 1/50. PubMed: 17103222WB: 1/1500. Detects a band of approximately 32 kDa (predicted molecular weight: 27 kDa).
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...