


Product Datasheet

Olfactory receptor 13C4 antibody (orb765893)


Description	Rabbit polyclonal antibody to Olfactory receptor 13C4
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OR13C4. AA range:201-250
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35
Expiration Date	12 months from date of receipt.



Western Blot
analysis of various
cells u...



Immunofluorescence
analysis of LOVO
cell...



Western blot
analysis of lysates
from CO...