



Product Datasheet

Olfactory receptor 2T2/35 antibody (orb765909)


Description	Rabbit polyclonal antibody to Olfactory receptor 2T2/35
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OR2T2/2T35. AA range:61-110
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
various
cells u...



Western
blot
analysis of
lysates
from HU...



Western
blot
analysis of
the lysates
fro...