


Product Datasheet

Olfactory receptor 4A15 antibody (orb765911)


Description	Rabbit polyclonal antibody to Olfactory receptor 4A15
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OR4A15. AA range:261-310
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36
Expiration Date	12 months from date of receipt.



Western Blot
analysis of various
cells u...



Immunofluorescence
analysis of A549
cell...



Western blot
analysis of lysates
from CO...