


Product Datasheet

BUB1 antibody (orb76592)

Description	Mouse monoclonal antibody to BUB1
Species/Host	Mouse
Reactivity	Human, Mouse, Other, Xenopus
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Fusion protein, corresponding to amino acids 281-419 of Human Bub 1.
Target	Bub1
Preservatives	Purified antibody (from supernatant) containing PBS 0.1% sodium azide;
Concentration	1.1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O08901 , O43683
Entrez	12235 , 296137 , 699
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try 1g/ml final concentration
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...