


Product Datasheet

Ca19.9 Sialyl Lewis a antibody (orb76593)

Description	Mouse monoclonal antibody to Ca19.9 Sialyl Lewis a
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, WB
Immunogen	Ten precipitin lines obtained after immunodiffusion using monoclonal antibody NS 19-9 and mucins isolated from an ovarian cyst of a OLe (a+b-) patient
Target	CA19-9
Preservatives	Purified antibody in phosphate buffered saline containing 10mM sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	Spleen cells from immunised BALB/c mice were fused with cells of the Sp2/0 myeloma line
Isotype	IgM
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try at 2ug/ml on FFPE tissue
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...