


Product Datasheet

Olfactory receptor 7C1 antibody (orb765947)


Description	Rabbit polyclonal antibody to Olfactory receptor 7C1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OR7C1. AA range:199-248
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
various cells
u...



Western blot analysis of lysates from He...



Western blot analysis of the lysates fro...