



Product Datasheet

p70 S6 kinase beta antibody (orb766014)


Description	Rabbit polyclonal antibody to p70 S6 kinase beta
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human p70 S6 Kinase beta. AA range:389-438
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42
Expiration Date	12 months from date of receipt.



Western Blot analysis of 293 cells using...



Western Blot analysis of various cells u...



Western blot analysis of lysates from 29...