


Product Datasheet

Cytokeratin antibody (orb76606)

Description	Mouse monoclonal antibody to Cytokeratin
Species/Host	Mouse
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	Rabbit esophageal keratins
Target	Cytokeratin 13
Preservatives	Purified antibody containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	P3 x 63 Ag8 myeloma cells;
Isotype	IgG1k
Clonality	Monoclonal
Uniprot ID	p13646, p08730
Entrez	3860, 16663
Dilution Range	Optimal antibody dilution should be determined by titration, however as guideline try at WB: 1/1000 - 1/3000. Predicted molecular weight: 50 kDa
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...