


Product Datasheet

Cytokeratin 3 antibody (orb76608)

Description	Mouse monoclonal antibody to Cytokeratin 3
Species/Host	Mouse
Reactivity	Bovine, Human, Rabbit
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	Rabbit corneal epithelial keratins
Target	Cytokeratin 3
Preservatives	Purified antibody containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	P3 x 63 Ag8 myeloma cells;
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P12035 , Q29426
Entrez	3850
Dilution Range	Optimal antibody dilution should be determined by titration, however as guideline try at WB: 1/1000 - 1/3000. Predicted molecular weight: 65 kDa.Flow Cyt: Use 2µg for 10 ⁶ cells.
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...