



Product Datasheet

PKC zeta antibody (orb766101)


Description	Rabbit polyclonal antibody to PKC zeta
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PKC zeta. AA range:526-575
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80
Expiration Date	12 months from date of receipt.



Western blot analysis of various cells u...



Western blot analysis of COS7 cells usin...



Immunofluorescence analysis of NIH/3T3 C...