


Product Datasheet

Desmoglein 2 antibody (orb76611)

Description	Mouse monoclonal antibody to Desmoglein 2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	N terminal amino acids 37-164 of Human Desmoglein 2 fused to a proprietary tag. Sequence used LLPKHPLVLRQKRAWITAPVALREGEDLSKKNPIAKIHSDLAEEERGLKITYKYTGK GITE PPFGIFVFNKDTGELNVTSILDREETPFLLTGYALDARGNNVEKPLELRIKVLDIN DNEPVFTQDVF
Target	Desmoglein 2
Preservatives	Purified antibody containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	NS1 myeloma
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q14126 , O55111
Entrez	13511 , 1829
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try at WB: 1/100. Detects a band of approximately 160 kDa (predicted molecular weight: 122 kDa).
Expiration Date	12 months from date of receipt.



Western blot analysis of HaCat keratinoc...